



UPPSALA
UNIVERSITET

What is learnt about SANS Instruments and Data Reduction from Round Robin Measurements? A Polymer Latex ‘Standard’

Adrian R. Rennie, Maja S. Helsing, Katy Wood,
Elliot P. Gilbert , Lionel Porcar, Ralf Schweins,
Charles D. Dewhurst, Peter Lindner,
Richard K. Heenan, Sarah Rogers, Paul D. Butler,
Jeff Krzywon, Ron E. Ghosh, Andrew J. Jackson




UPPSALA
UNIVERSITET

What is standardisation?

Instrument	AAA
User	BBB
Software	CCC
Date	DDD
Sample	SSSS

Instrument	WWW
User	XXX
Software	YYY
Date	ZZZ
Sample	SSSS



Do I get the same result? Has the sample changed?
How sure am I?

How do we obtain similar results?



UPPSALA
UNIVERSITET

Are results consistent?

- Is the size (distribution) the same as that from electron microscopy, light scattering, GPC ?
- Does SAXS and SANS give the same result?
- Do I have the same conclusion from model fitting and inversion procedures?


Do we understand the differences?



UPPSALA
UNIVERSITET

More than Calibration

- Wavelength
- Distance
- Angle
- Intensity
- Resolution
- Uniformity of detector
- etc.



How do I check these quickly?



UPPSALA
UNIVERSITET

Different Questions?



User: Do I understand the data?
Are my results publishable?



Instrument scientist:
Why are results different? Can the user publish the data?



**Facility
Manager:**

My instruments are the best?

Everyone needs to understand better!



UPPSALA
UNIVERSITET


Why Standardisation?

Comparisons:

- Samples
- Instruments
- Procedures
- Techniques
- Software

**Provide understanding of
small-angle scattering!**


**Co-operation and comparison
helps this understanding**





UPPSALA
UNIVERSITET

Instruments




NG7, NCNR, USA



Bragg Institute, Australia



SANS2D –ISIS, UK



D22 and D11, ILL, Grenoble, France






UPPSALA
UNIVERSITET


Round Robin Sample

PS3 Polystyrene latex in D₂O



SEM 5 keV uncoated
latex on Si wafer

8% - Diffracts light



Static light scattering – ALV
HeNe laser $R_g = 56 \text{ nm}$

$R = 716 \text{ \AA} \pm 2 \text{ \AA}$




UPPSALA
UNIVERSITET

Orders on
surfaces and in
the bulk

Topographic
mode
Nanoscope SM


AFM






UPPSALA
UNIVERSITET

Dynamic Light Scattering



$$R_H = 71 \text{ to } 73 \text{ nm}$$

Polydispersity $\sim 5\%$




Good agreement with
SLS and SEM

$\lambda = 633 \text{ nm}$, $\theta = 90^\circ$: data measured also at other angles



UPPSALA
UNIVERSITET

Scattering from Spheres



Form Factor: $Q / \text{\AA}^{-1}$

$$P(QR) = \{3(\sin(QR) - QR \cos(QR))/(QR)^3\}^2$$

Intensity:


$$I(Q=0) = (\Delta\rho)^2 n V^2$$



UPPSALA
UNIVERSITET

0.43% Latex in D₂O

1 mM NaCl






UPPSALA
UNIVERSITET

Differences – Measured data

0.43% Latex in D₂O 1mM NaCl




Are some
data wrong?



UPPSALA
UNIVERSITET

Differences – Measured data

0.43% Latex in D₂O 1 mM NaCl



Are some
data wrong?



UPPSALA
UNIVERSITET

Presenting Data


Conclusions

- Logarithmic scales do not show everything well!
- Data are not necessarily wrong but perhaps misinterpreted
- Need more information – better description of metadata and uncertainties





Simple Fits – SasView Spheres


UPPSALA
UNIVERSITET



Polydispersity: 8%



3%




SANS2D data: Which fit is better? Both show systematic deviations!



UPPSALA
UNIVERSITET

Which fit is best?

- Better (when choosing from 2) but neither is best!
- 8% polydispersity has smaller χ^2 but misses all large Q features
- Need more information



Fit with 8% polydispersity

R either 687 Å or 703 Å

(polydispersity 8% or 3%)



UPPSALA
UNIVERSITET

Model Fitting

Need to include:

Resolution

Polydispersity

Multiple scattering

Interactions ?

Effects are similar but not identical

Variation with Q and concentration is different




UPPSALA
UNIVERSITET

Different Concentrations

$$R = 721 \text{ \AA}$$

$$\sigma = 21 \text{ \AA} (\sim 3\%)$$



Sample	Concentration / wt %
A	0.43
B	0.11
C	0.03


D22 data: simultaneous fit hs2m

includes resolution and double scattering



UPPSALA
UNIVERSITET

Monte Carlo Simulation




D22 data MC simulated with NCNR IGOR programs (J. G. Barker, S. G. Kline et al)



UPPSALA
UNIVERSITET

Multiple Scattering




Data: D22, ILL 12 Å



UPPSALA
UNIVERSITET

Compare Ratio - Data & MC



Monte Carlo modelling can account for smearing by multiple scattering
Calculations for R = 705 Å 4% polydispersity



UPPSALA
UNIVERSITET


Analysis Methods

- Guinier analysis – limited fit information and needs low Q – no resolution
- Modelling scattering – multiple data sets and detailed knowledge of instrument/resolution needed. Only limited multiple scattering.
- Monte Carlo – needs precise instrument geometry. Background is difficult but MC can include coherent multiple scattering



UPPSALA
UNIVERSITET

Inversion to $g(r)$





Inversion works when data extends over adequate Q range but (a) resolution is rarely considered and (b) $g(r)$ needs to be interpreted.



UPPSALA
UNIVERSITET

Challenges


Missing/incorrect
metadata



Resolution not Gaussian



$\delta\lambda$ changes
with $\delta\theta$



Detector normalization
changes with configuration



UPPSALA
UNIVERSITET

Challenges



- Incorrect or missing metadata (e.g. to calculate resolution or multiple scattering)
- Wavelength resolution with velocity selector can depend on collimation
- Resolution function may not be a Gaussian – particularly on ToF instruments
- Detector normalization: inappropriate ‘flat field’ can distort data



Conclusions – What have we learnt?


UPPSALA
UNIVERSITET

Data that can
be modelled
reliably helps
comparisons



Compare instruments and software

Systematic deviations
are often the largest
source of uncertainty in
interpretation



ToF and const λ measurements
provide beneficial comparisons



UPPSALA
UNIVERSITET

Conclusions

- Well-known form factor has identified problems with resolution, detector normalization and software
- Single wavelength data were easier to model in detail
- Time-of-flight SANS data with a wide Q-range and with good resolution highlighted multiple scattering
- **Systematic deviations are often the largest source of uncertainty in interpretation**
- Many other 'unknown' samples can show similar effects



UPPSALA
UNIVERSITET

Recommendations

- Regular comparisons of instruments and procedures as well as software are helpful
- Data formats and publishing standards need to include uncertainty from systematic effects as well as counting statistics
- Do not be tempted to scale data to ‘match’ without allowing for resolution!
- Descriptions of data are essential - e.g. how is resolution described, σ , FWHM etc.?
- Density matched ‘sealed’ sample for long term reproducibility would be helpful



UPPSALA
UNIVERSITET

Thanks

- Facilities and the Funding Agencies for the facilities
- Co-operation between many instrument scientists
- www.cansas.org

Thank you for listening



UPPSALA
UNIVERSITET

Contact

Join in these activities?

Adrian.Rennie@physics.uu.se

Lunch time canSAS session at 12.40 today

